



All Blessings Flow
EQUIPMENT & SUPPLIES CATALOG

Please call 434-422-8888 to see if items are in stock.



Contents

- 3 Dressings Aids
- 5 Specialty Clothing
- 6 Personal Care
- 8 Bathing Accessories
- 10 Toileting
- 12 Dining Accessories
- 13 Beds and Accessories
- 16 Visual Aids
- 17 Household
- 18 Transfer Assistance
- 21 Wheelchairs
- 22 Orthopedic
- 24 Exercise Equipment
- 26 First Aid
- 27 Health Monitoring
- 28 Breathing Equipment and Accessories
- 29 Communication Devices
- 29 Safety Devices
- 30 Medication Accessories
- 30 Fine Motor Skills Tool



Dressing Aids

BUTTON HOOK

A handheld wand with a wire loop at the end used to pull a button through a buttonhole using only one hand.

DRESSING STICK

A long wooden or plastic stick with a hook on one end that can assist in various dressing tasks such as grabbing a shirt or pants.

LEG LIFTER

A long piece of stiff fabric with a looped end used to move a leg from one area to another.

SHOEHORN – LONG HANDLE

A contoured piece of hard plastic used by placing it in the heel of a shoe, allowing an individual to slip their foot into the shoe with minimal effort.



SOCK AID

Assists an individual in putting on socks by providing an open space for the foot to go through, allowing an individual to pull the sock onto the foot and over the heel in with ease.



REACHER

A long plastic or metal stick with a clamp at the end that is used to pick up items when the individual is lacking range of motion, cannot bend, or reach forward.

Specialty Clothing



COMPRESSION SOCKS

A special type of hosiery made to help in the prevention of swelling and venous clotting in the lower legs. Various heights are available including below knee or above knee.



HOSPITAL SOCKS

A sock with non-slip tread on the bottom that helps prevent falls and a non-constrictive top that won't be tight around the leg.



HOSPITAL GOWN

Short, collarless gown that ties in the back which allows for easy application and removal of clothing.

Personal Care



ADULT PULL-ON UNDERWEAR

Disposable adult underwear that offers protection for incontinence.



BRIEFS WITH TABS

Offers protection for incontinence and can easily be taken on and off utilizing tabs.



BELTED SHIELDS

An incontinence pad with reusable straps that help keep it in place during activities.



INCONTINENCE PADS

A pad placed in the underwear that offers protection for incontinence throughout the day.



SHAMPOO CAP

A ready-to-use cap that is placed on an individual's head and then massaged to shampoo and condition hair without any rinsing required.



SKIN CREAM PROTECTANTS

Used to prevent and treat rashes, burns, scrapes, and many other skin irritations by forming a protective layer on the skin.



WIPEs

Moistened wipes used to quickly clean an individual's skin and sensitive areas.



BATH WIPEs

Single-use, moistened clothes that cleanse and moisturize skin without any need for rinsing. Useful for individuals who are unable to bathe.

Bathing Accessories



SHAMPOO BASIN - INFLATABLE

Allows for portable, in-bed shampoo washing and cleansing. The inflatable basin has a secure frame made for the head and shoulders to rest in.



SHAMPOO TRAY

Provides a way for comfortable hair washing from a seated position with assistance from a caregiver. Stability straps may be available with some items.



WASH BASIN

A portable plastic bin used for cleaning person or equipment.



SHOWER CHAIR

A chair with rubberized feet and a durable, plastic seat used to provide additional support in the bath or shower. Some models also have back support.



SHOWER CHAIR – TILT-IN-SPACE

A rolling shower chair for easy access to bathing. It can be rolled out of the bathroom without the need to transfer the person on a wet surface.



TUB CHAIR LIFT - ELECTRIC

A device used to help individuals with mobility deficits transfer into and out of the bathtub. Electric power relieves the caregiver of having to provide strenuous physical assistance.



TUB SAFETY RAIL

A rail that clamps onto the edge of the bathtub and is securely tightened to help while transferring into and out of the tub.



TUB TRANSFER BENCH

An extra-long shower bench that is placed over the tub wall, allowing an individual to sit down in the chair while outside of the tub and then slide across the chair in order to transfer into the tub.

Toileting



BEDSIDE COMMODE

Adult chair used for toileting that is made of a frame equipped with a toilet seat and a removable bucket.



BEDSIDE COMMODE - DROP ARM

A bedside commode with easy-release drop arms for more convenient side transfer in tight spaces.



RAISED TOILET SEAT

Allows individuals to sit up higher while using the toilet and makes it easier for them to return to standing when they are finished.



TOILET SAFETY RAILS

A frame that attaches to the toilet to provide support while sitting down and getting up. Height is adjustable to accommodate for a raised toilet seat.



SELF-WIPE TOILETING AID

An extended plastic wand with a small clip at the end used to hold toilet paper. Once toilet paper is secured on the end of the wand, it may be used during toileting tasks.



BEDPAN

A receptacle used for toileting of bed-bound individuals. Easy to clean and allows for relief without having to leave the bed.



URINAL - HANDHELD

A portable container for males to use while traveling or unable to get up to use the facilities to urinate.

Dinning Accessories



ADAPTIVE TUBING

Foam tubes used with dinning utensils, pens, and toothbrushes to provide a wide, secure grip for those with unsteady hands.



BIB

A protective covering over the chest area used during eating, especially for individuals with decreased upper body strength.



UTENSILS WITH BUILT-UP HANDLES

Specialized utensils with a wider grip that allow individuals to eat more efficiently if they have trouble using their arms or hands.

Beds and Accessories



BED RAIL

An adaptation installed on a bed to prevent an individual from falling out or to grasp onto for repositioning.



BED TRAPEZE

Attaches to the bed to function as a grab bar, providing leverage for individuals so they can position themselves without help.



BEDSIDE TABLE – HOSPITAL STYLE

An adjustable, wheeled table that can easily be moved around and over the bed.



HOSPITAL BED

A fully electric hospital bed that easily allows for adjustments to the head and foot of bed as well as the height of the frame with the push of a button.



AIR MATTRESS SYSTEM

A mattress with air cells and sensors that ensure weight is well distributed for maximum comfort and prevention of pressure injuries.



ALTERNATING PRESSURE MATTRESS PAD

An inflatable mattress topper with a quiet pump that automatically provides variable waves of pressure to prevent pressure injuries.




FOOT TENT – BLANKET SUPPORT

A structure that attaches to the end of the bed in order to prevent blankets from creating too much pressure on an individual's foot and allowing air flow for the foot. Commonly used for individuals with sensory issues or injuries to the area.



HEEL PROTECTOR

Used by individuals who spend most of their time lying down to relieve pressure, prevent twisting of the leg, and avoid friction to prevent pressure ulcers.



POSITIONING WEDGE

A sturdy foam wedge used for support and stabilization that is useful when re-positioning an individual to prevent pressure injuries.

Visual Aids



MAGNIFIER – FLOOR STAND

A large magnification lens attached to a movable arm that allows an individual to see reading materials or fine hand work.



MAGNIFIER – HANDHELD

A portable lens through which the user may view small objects or text.



READING GLASSES

Glasses with magnifying lenses that aid the wearer in seeing text or small objects more closely.

Household



GERI CHAIR

A large, padded reclining chair with wheels designed to allow patients recovering from illness and surgery to get out of bed and sit comfortably.



STEP STOOL

Used to assist in accessing out-of-reach objects such as high shelves. Can also be used for step-on/step-off exercise to improve strength.



BED TRAY

A tray with foldable legs that fits over one's lap and functions as a small and easily portable table.



LAUNDRY CART

A metal cart on wheels that can be used for transporting laundry, groceries, or other heavy loads. Is collapsible for easy storage.

Transfer Assistance



SIT-TO-STAND LIFT

A sit-to-stand manual lift with hydraulic pump that provides transfer support that is ideal for individuals who are partially weightbearing. Wheels and low clearance design allow for easy maneuverability around the home and under furniture.



HOYER LIFT – MANUAL OR ELECTRIC

A lift used to safely move an individual from one surface to another. The manual style is operated by a handle and hydraulic pump that moves the sling-bearing arm up and down. The electric style is equipped with a motor that powers the arm movement.



STAIR LIFT

A device designed for moving an individual up and down the stairs. A chair or platform is attached to a mounted rail that a person then rides for safe and easy ascent or descent.



SECURITY POLE – GRAB BAR

A floor to ceiling pole with a pivoting curved bar used for balance and ease of sitting and standing. Easily installs via tension without extra tools.



PREVALON TURN AND POSITION SYSTEM

Provides a safe way to move and rotate individuals in bed. The system stays beneath the person which makes it always ready and easy to use.

**EASY UP SEAT ASSIST**

A device that gives individuals a boost when going from seated to standing. No electricity needed.

**GAIT BELT**

A safety device used to support individuals as they move and prevent falls. It also decreases a caregiver's risk of back injury while helping a person and provides a safe, stable surface to hold onto during movements.

**SWIVEL SEAT CUSHION**

A soft, flexible cushion that is placed on a seat and swivels to make getting in or out of the car, wheelchair, or chair easier and more comfortable.

**TRANSFER BOARD**

A piece of equipment that helps transfer an individual from one surface to another as they slide across.

Wheelchairs



STANDARD WHEELCHAIR

A manual wheeled mobility device in which the user sits and can propel themselves by turning the large wheels by hand.



TILT-IN-SPACE WHEELCHAIR

A manual wheelchair with adjustable armrests, backrest, and push-bars that both reclines and tilts to accommodate the particular needs of an individual.



RECLINING WHEELCHAIR

A manual wheelchair with a seatback that may be adjusted to many different angles and padded arm rests that may be removed or extended.



TRANSPORT WHEELCHAIR

A lightweight wheelchair that is smaller and easier to move around and maneuver through smaller spaces than a standard wheelchair. It is intended to be pushed by a caregiver as the small wheels don't permit the user to propel him/herself.

Orthopedic



21

ANKLE BRACE

Provides support for the injured or unstable ankle. Is slim enough to fit into a shoe for use during activities.



KNEE IMMOBILIZER

A brace that prevents the knee from any movement by keeping the knee in one position to support healing after an injury.



22

NIGHT SPLINT

A foot splint with easy-on straps that is used for the treatment and relief of pain associated with plantar fasciitis by providing a stretch to the bottom of the foot.



WALKING BOOT

A boot that restricts movement in the lower leg and ankle to facilitate healing after an injury or surgery. It may be worn after surgery to reduce risk of reinjury.



SHOULDER BRACE

Prevents excess movement to aid in recovery after an injury. Also supports the shoulder by promoting appropriate posture.



WRIST BRACE

Maintains the position of an individual's wrist and facilitates healing for the area. It inhibits movement that causes pain or discomfort and may also prevent further injuries.



POLAR CARE CUBE

A cold therapy system used to reduce pain and swelling in the shoulder or knee. A cooler is filled with ice and water that flows through a slim tubing system to be spread out around the targeted area.



TENS UNIT

A transcutaneous electrical nerve stimulation (TENS) unit uses two pads that send a small electrical current to a specific area of the body to provide pain relief.

Exercise Equipment



ANKLE WEIGHTS

Light weights with straps that can be added to walking or targeted leg exercises to increase strength and endurance.



EXERCISE BIKE

A stationary bicycle used for strengthening of the thigh and calf muscles and facilitating range of motion.



EXERCISE MAT

A foam mat used as a base to perform exercises and stretches on. It can increase comfort while preventing injury when performing exercises on the floor.



FOAM ROLLER

Used by rolling the device over different muscle groups to relieve soreness and muscle tightness.



THERABAND

Latex resistance bands that are used for strength training or home exercise programs to facilitate healing, stretching, and strengthening.



THERAPUTTY

Special putty that is used to provide light resistance to increase grip strength, fine motor skills, dexterity, and increased finger and hand strength.

First Aid



ACE WRAP

A flexible elastic bandage that provides compression and support. Different sizes available.



DRESSINGS - STERILE

4x4 or 2x2 gauze squares in individual packaging used for protection or cleansing of a wound.



HEATING PAD

A square pad that plugs into an outlet to provide adjustable heat therapy.



HOT/COLD PACKS

Small packs that are either frozen or heated then used to provide therapeutic relief.

Health Monitoring



BLOOD PRESSURE MONITOR

A digital machine that allows the user to check blood pressure with an automated cuff and the push of one button.



BLOOD PRESSURE CUFF AND STETHOSCOPE

Allows a caregiver to manually measure an individual's blood pressure.



GLUCOSE MONITOR

A portable small machine that measures the amount of glucose (blood sugar) in the blood. Especially used by people with diabetes.



PULSE OXIMETER


A small device placed on the finger that measures the heart rate and amount of oxygen in the blood.

Breathing Equipment and Accessories



CPAP MACHINE

An effective treatment for moderate to severe sleep apnea. It attaches to the user's face and nostrils, allowing for assistive breathing during sleeping hours.



NEBULIZER MACHINE

A machine that turns liquid medication into a mist so that it may be easily breathed into the user's lungs through a mask or mouthpiece.



OXYGEN TANK CART

A small cart with wheels that allows for the easy transportation of an oxygen tank. *Oxygen tank not included.*

Communication Devices



23

© Copyright 2010 Tellex, Inc.

ABLEPHONES – 7000VC

A landline phone with an easy to understand layout including large numbers and symbols, a lightweight corded phone, and a voice-activated speakerphone for users with auditory impairments and/or hand weakness.



HEARING AID

In-ear auditory support for individuals with hearing impairments. A microphone picks up sound and amplifies it for the user.

Safety Devices



ALARM SENSOR

A device that is attached from an individual's clothing to a bed or chair that sounds an alarm if the individual arises.



24

LIFE ALERT

An emergency button that one carries so when an emergency situation such as a fall or severe accident occurs, the user may press the button and emergency services will be sent to the residence.


Medication Accessories



PILL CRUSHER

Used to crush tablet medications into a powder form that can be mixed into another substance for easy swallowing.

**Check with your medical provider or pharmacist before crushing prescribed medications.*



WEEKLY PILL ORGANIZER

Allows the user to divide medications per day and time (morning/evening).

Fine Motor Skills Tool



25

ACTIVITY APRON

An apron with many pockets, specific fabrics, and items attached. Allows for practice of sensory skills as well as maintenance of hand motion skills

Please note:
Items in photos
are for example
only. Style of
exact items
available may
differ.



<https://www.allblessingsflow.org/>

Image References:

1. <https://www.silverts.com/show.php/product/19140-womens-mens-hospital-socks-non-skid-resistant-anti-slip-resistant-grip-socks-fuzzy-gripper-socks-regular-size-xl-bariatric/grey/alt-web>
2. <https://www.carewayclub.com/shop/incontinence/incontinence-briefs-underwear/mckesson-regular-absorbent-underwear>
3. <https://www.carewayclub.com/shop/incontinence/incontinence-belted-and-knit-products/prevail-belted-shields-incontinent-undergarment>
4. <https://cleanlifeproducts.com/hair-care/>
5. <https://www.cascadehealthcaresolutions.com/medline-readybath-total-body-cleansing-washcloths-p/msc095305.htm>
6. <https://www.cmecorp.com/fei-inflatable-shampoo-basin.html>
7. <http://www.accessiblehomesadvisor.com/bathtub-lifts-for-disabled/>
8. <https://www.amazon.com/Drive-Medical-Transfer-Adjustable-Backrest/dp/B002VWK0T6>
9. <https://www.amazon.com/Medlines-Guardian-Toilet-Adjustable-Bathroom/>
10. <https://thetoiletzone.com/best-toilet-tissue-aid/>
11. <https://www.bluechipmedical.com/mattress-systems/air-mattress/supreme-air/>
12. https://www.amazon.com/dp/B081B7SZ57/ref=cm_sw_em_r_mt_dp_U_2dhLEbCCVYSZK
13. <https://www.neppt.com/products/adjustable-blanket-support-lifters-bed-foot-tent-blankets-raiser-heavy-duty-steel-medical-accessories>
14. <https://www.alimed.com/alimed-body-positioning-wedges-uncovered>
15. <https://medmartonline.com/patient-room/geri-chairs/>
16. <https://www.amazon.com/Winsome-Wood-98821-Stockton-Natural/dp/B001U88698>
17. <https://www.amazon.com/Whitmor-Deluxe-Utility-Extra-Large/>
18. http://www.invacare.com/cgi-bin/imhqprd/inv_catalog/prod_cat_detail.jsp?s=0&prodID=GHS350&catOID=-536885361
19. <https://stander.com/product/1100-security-pole-curve-grab-bar/>
20. <https://carex.com/products/carex-upeasy-seat-assist-plus>
21. <https://www.mountainside-medical.com/products/procare-lace-up-ankle-brace>
22. <https://www.medicaldepartmentstore.com/Dorsal-Night-Splint-Plantar-Fasciitis-Foot-Drop-p/66113.htm>
23. <https://teltex.com/ap-7000-voice-dialer/>
24. <http://www.lifealert.com/>
25. https://www.alimed.com/discovery-apron.html?pid=97441&gclid=CjwKCAjwpqv0BRABEiwA-TySwSjdgigp2wzuYpOkNT0Radh2kARCuWaYzClZPRIfagqrh6MmtV4XghoCCwwQAvD_BwE